

# Unbreakable™

*Unbeatable strength*



- › Nylon thermoplastic
- › Aesthetic
- › **Guaranteed unbreakable**
- › Lightweight yet durable
- › Ideal degree of flexibility
- › Stain resistant
- › Repairable and rebasable
- › Biocompatible
- › Millions of successful cases prescribed





## BENEFITS

Ideal for partial dentures and unilaterals  
 Economical, aesthetic solution for missing teeth  
 Metal-free  
 Gingival veneers to mask exposed root surfaces  
 Ideal for patients allergic to monomer or nickel  
 Used as a transitional restoration during the healing period for implants  
 Presence of tori, oral cancer, cleft palate



## TIPS

1. For best results use lightweight alginate impression materials.
2. Use TCS Green Points for making small adjustments and TCS Brown Rubber Points to smooth adjusted areas.
3. For a touch-up polish, Green Hi-Shine polishing paste (item #4209-02) is recommended.
4. Soak the partial in very hot clean tap water for about 1 minute before initial insertion. When case cools to a tolerable temperature gently insert it into patient's mouth. (This is only recommended for 1<sup>st</sup> insertion)
5. Recommend TCS Dental Appliance Cleaner to your patients for optimal maintenance.
6. Use the TCS Animation Video (1.5 min.) & Health Briefs with Terry Bradshaw (5 min.) for patient chairside education.  
 Visit [www.tcsdentalinc.com](http://www.tcsdentalinc.com) to download videos.



## SHADES



Light Pink | Standard Pink | Light/Dark Pink | Dark Pink | Natural (clasps)



## INSURANCE CODES

D5225 Maxillary Partial Denture  
 D5526 Mandibular Partial Denture

*"I love my TCS Partial Denture. It is comfortable, easy to care for, and best of all no one can tell I am wearing it."*

*~Violette, 72 years old*

Specify Unbreakable™ on your Lab Slip

