

Karadent™ by tcs®

The Modern Fit



For upgraded full dentures, semi-flexible partials, clear clasps, frameworks & nightguards

- › Semi-rigid, microcrystalline polymer thermoplastic
- › No volume shrinkage for an exceptional fit
- › Minimal water absorption
- › Easy to repair, rebase, and reline with acrylic
- › An aesthetic upgrade from acrylic dentures
- › Available in shade Clear for complete transparency
- › A modern fit, lightweight and durable



Monomer Free
Hypoallergenic



INDICATIONS & BENEFITS

Ideal for full dentures, clear clasps, frameworks, nightguards & partial dentures when more stability is needed
 Economical, aesthetic solution for missing teeth
 Metal-free solution
 Ideal for patients with allergies to monomer and/or nickel
 Transitional restoration during the healing period for implants



SHADES



Standard Pink



Dark Pink



Clear



INSURANCE CODES

D5225 Maxillary Partial Denture
 D5526 Mandibular Partial Denture
 D5110 Maxillary Complete Denture
 D5120 Mandibular Complete Denture



TIPS

1. For best results use lightweight alginate impression materials.
2. Karadent™ is quick and easy to adjust. It can be adjusted using traditional acrylic tools and techniques.
3. We recommend using the TCS Pink Zirconia Bur (5035Z) or any carbide bur for minor adjustments and the TCS Brown Rubber Points (4513-01, 4514-01) to smooth the adjusted areas. Use light pressure in a scraping motion.
4. For a quick re-polish apply Green Hi-Shine polishing paste (4209-02) with a dry Mini Rag Buff (4221-01).
5. **Recommend the tcs® Concentrated Dental Appliance Cleaner to your patients for optimal appliance maintenance.**

Chairside adjustment tools and FREE promo material available at tcsdentalinc.com.

Karadent's unique formula achieves the unbreakable strength of nylons with the simplicity of acrylics to create the most stable and durable removable appliances the industry has been waiting for.

Specify Karadent™ on your Lab Slip

